

DILUENT PH7



Wash solution for erythrocytes

Diluent PH7 is an isotonic phosphate buffered saline solution with pH 7,0 which is produced for use at blood centers.

To extend the shelf life a preservative is added.

Diluent PH7 is provided in packaging of 1 liter, 5 liter and 10 liter.

Diluent PH7 10 liter

TV001

Diluent PH7 5 liter

TV002

Diluent PH7 12 x 1 liter

TV003-12


ISO 9001:2015 certified



labexreagents.com

LABEX

OF SCANDINAVIA



For further information, please contact us:
LABEX REAGENS AB • Ekslingan 6 • SE-254 67 • Helsingborg • Sweden
phone:+46 (0)42 32 40 00 • e-mail: reagents@labex.com web:labexreagens.com

LABEX was founded 1973 and is today a successful supplier of diagnostics and instruments for laboratories in hospitals, universities and pharmaceutical industry.

LABEX has produced cellcounting reagents for more than 20 years and has a well established world wide network of distributors.

LABEX is ISO 9001:2015 certified and all products are CE-marked to the IVD directive.



LE20010 02/2020



labexreagens.com

LABEX

OF SCANDINAVIA